

Transcranial Extraocular Light

**** WARNING : NEVER LOOK DIRECTLY INTO THE LIGHT ****

CAUTION: For use with any contact51 Frequency Generator only

Usage:

1. Plug into the **PROBES** socket on the generator
2. Insert earpieces one into each ear
3. Turn generator on
4. Use for any length of time as required

SPECIFICATION:

Size: 31(H) x 12(w) x 8(D) millimetres per earpiece

Ear tube diameter:10(D) x 20(L)

Weight: ~ 22g

- Power: Contact51 Frequency Generator **only** (Probes socket)

Leads: 500mm piece-piece; 1.0M connecting with 3.5mm jack plug

Output colours available:

Red ~700nm

Green ~530nm

Blue ~510nm

Amber ~620nm

White covers all above - visible colour spectrum

We manufacture equipment and are not practitioners. We are therefore unable to provide any advice. Before using any equipment for medical related issues you should always first consult your medical professional.